

## TURNAROUND TIMES FOR STAT CRITICAL LAB TESTS

The following menu of Laboratory tests may be ordered STAT at all times if rapid communication of the results is required. The stated turnaround time in the laboratory is defined as the time interval from specimen **RECEIPT** in the lab to result entry.

STAT Chemistry Panels/Screens		
Test Code	Test Description	Turn-Around Times (minutes)
NA, K, CL	Electrolytes	40
BMP	Basic Metabolic Panel	40
CMP	Complete Metabolic Panel	40
LIPPR	Lipid Profile	40
HEPFP	Hepatic Function Panel	40
RENAL	Renal Panel	40
DAU	Drugs Abuse Screen	60
CKCKMB, TROPN, PBNP	Cardiac Screens	50
TSH, FT4	Thyroid Screens	60

Other STAT Chemistry Tests		
Test Code	Test Description	Turn-Around Times (minutes)
AMMON	Ammonia	40
AMYL	Amylase	40
BHCG	Beta-Human chorionic gonadotropin	75
BHB	Beta-Hydroxybutyrate	40
TBIL DBIL	Total/Direct Bilirubin	40
CK	Creatine Kinase	40
CRP	C-Reactive Protein	40
CSFGTP	CSF Glucose/Total Protein	40
IONCA	Ionized Calcium	40
LACT	Lactic Acid	40
LIPA	Lipase	40
MG	Magnesium	40
OSMOLS (Serum) OSMOLU (Urine)	Osmolality	40
PHOS	Phosphorus	40
PREGS (Serum) PREGU (Urine)	Pregnancy	40
URIC	Uric Acid	40

STAT Drug Tests		
Test Code	Test Description	Turn-Around Times (minutes)
ACETA	Acetaminophen	40
ALCO	Alcohol	40
DIG	Digoxin	40
SALI	Salicylate	40
VANC	Vancomycin	40

STAT Coagulation Tests		
Test Code	Test Description	Turn-Around Times (minutes)
PROI	Prothrombin Time	40
PTT	Partial Thromboplastin Time	40
DDIMR	D-Dimer	60
FBG	Fibrinogen	60
*RAPTEG	Rapid Thromboelastogram	75

STAT Hematology Tests		
Test Code	Test Description	Turn-Around Times (minutes)
HH	Hemoglobin/Hematocrit	15
CBC (with no diff)	Complete Blood Count with no differential	15
CBCD	Complete Blood Count with differential	40
CCT/FCT	CSF/Fluid Cell Count	60
RET	Reticulocyte Count	20
UAC	Urinalysis with Culture if indicated	30

STAT Microbiology Tests		
Test Code	Test Description	Turn-Around Times (minutes)
GS, (A)ACGS, (A)AGWND, AGSPUT	Gram Stain (Only Gram Stain result can be STAT)	45
WP	Wet Prep	30
*MEP	Rapid Meningitis/Encephalitis Panel	180
RSV	Respiratory Syncytial Virus Screen	45

STAT Blood Bank		
Test Code	Test Description	Turn-Around Times (minutes)
XM	Crossmatch	60

For STAT testing requested from physician's offices, please allow a minimum of **one hour** for transport of specimens within the medical district to the hospital laboratory.

Any test **ONLY** offered at HH Main Laboratory is indicated with \*. Please allow a minimum of **one hour** for transportation of specimens to the HH Main Laboratory from Madison Hospital.